

Anti-Rig-I/Ddx58 Monoclonal Antibody For Wb, Elisa - 095786

Item Number: CM0009681



Introduction

Anti-Rig-I/DDX58 rabbit monoclonal antibody (SWISS 095786) validated for WB and ELISA. Detects human RIG-I, a cytoplasmic RNA sensor pivotal in innate antiviral immunity. Suitable for infectious disease and interferon signaling research studies.

[Learn More](#)

Parameter	Value
Product Name	Rig-I/DDX58 Rabbit mAb
Remarks/Alias	RIG1; DDX58; RIG-I; RLR-1; SGMRT2; Rig-I/DDX58
Species	Human
Gene ID	23586 (Human)
Immunogen	Synthetic Peptide: A synthesized peptide derived from human Rig-I/DDX58
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	095786
Protein Weight	107kDa
Shipping	Ice bag