

Anti-Rig-I/Ddx58 Rabbit Monoclonal Antibody For Wb, Elisa - O95786

Item Number: CM0007108



Introduction

Rabbit monoclonal antibody targeting human Rig-I/DDX58 (RIG-I, DDX58). Validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, rat. Suitable for innate immunity and antiviral response research.

[Learn More](#)

Parameter	Value
Product Name	Rig-I/DDX58 Rabbit mAb
Remarks/Alias	RIG1; DDX58; RIG-I; RLR-1; SGMRT2; Rig-I/DDX58
Species	Human
GeneID (Human)	23586
GeneID	23586
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95786
Protein Weight	107 kDa
Shipping	Ice bag