

Anti-Rig-I/Ddx58 Rabbit Polyclonal Antibody For Wb And Elisa - O95786

Item Number: CM0014568



Introduction

CamelBio offers a specific rabbit polyclonal antibody against human Rig-I/DDX58, validated for Western blot and ELISA. Detects human and rat Rig-I, a crucial innate immune receptor for viral RNA sensing and antiviral signaling.

[Learn More](#)

Parameter	Value
Product Name	Rig-I/DDX58 Rabbit pAb
Remarks/Alias	RIG1; DDX58; RIG-I; RLR-1; SGMRT2; Rig-I/DDX58
Species	Human
GeneID	23586
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 825-925 of human Rig-I/DDX58 (NP_055129.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O95786
Protein Weight	107 kDa
Shipping	Ice bag