


Anti-Rims2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uq26

Item Number: CM0017932



Introduction

High-quality rabbit polyclonal antibody against RIMS2, suitable for WB and ELISA applications with cross-reactivity to human and rat. Ideal for studying synaptic exocytosis and related neurological pathways.

[Learn More](#)

Parameter	Value
Product Name	RIMS2 Rabbit pAb
Remarks/Alias	OBOE; RIM2; CRSDS; RAB3IP3; RIMS2
Species	Human
Gene ID (Human)	9699
Gene ID	9699
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human RIMS2 (NP_055492.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9UQ26
Protein Weight	160kDa
Shipping	Ice bag