

Anti-Rin2 Polyclonal Antibody For Wb And Elisa - Q8Wyp3

Item Number: CM0013462

Introduction

Rabbit polyclonal antibody targeting human RIN2 (Ras and Rab interactor 2), also known as RASSF4. Validated for WB and ELISA, cross-reacts with mouse. Ideal for endocytosis pathway studies.

[Learn More](#)

Parameter	Value
Product Name	RIN2 Rabbit pAb
Remarks/Alias	MACS; RASSF4; RIN2
Species	Human
Gene ID	54453
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8WYP3
Protein Weight	100kDa
Shipping	Ice bag