


## Anti-Ring1B/Rnf2 Monoclonal Antibody For Wb, Ip, And Elisa - Q99496

Item Number: CM0009853



### Introduction

High-quality rabbit monoclonal antibody targeting human RING1B/RNF2 (RING2). Validated for WB, IP, and ELISA, with cross-reactivity to mouse and rat. Essential for studying Polycomb repression and histone H2A ubiquitination.

[Learn More](#)

| Parameter        | Value                                                        |
|------------------|--------------------------------------------------------------|
| Product Name     | RING1B/RNF2 Rabbit mAb                                       |
| Remarks/Alias    | BAP1, DING, BAP-1, HIP13, RING2, LUSYAM, RING1B, RING1B/RNF2 |
| Species          | Human                                                        |
| GeneID (Human)   | 6045                                                         |
| GeneID           | 6045                                                         |
| Immunogen        | Synthetic peptide (aa 100-200 of human RING1B/RNF2, Q99496)  |
| Source           | Rabbit                                                       |
| Category         | Monoclonal Antibodies                                        |
| Application      | WB, IP, ELISA                                                |
| Cross Reactivity | Human, Mouse, Rat                                            |
| SWISS            | Q99496                                                       |
| Protein Weight   | 38kDa                                                        |
| Shipping         | Ice bag                                                      |