


Anti-Ring1B/Rnf2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa, Chip - Q99496

Item Number: CM0030750



Introduction

Rabbit polyclonal antibody against RING1B/RNF2 (E3 ubiquitin-protein ligase RING2), validated for WB, IHC-P, IF/ICC, IP, ELISA, ChIP. Cross-reacts with human, mouse, rat. Ideal for chromatin and epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	RING1B/RNF2 Rabbit pAb
Remarks/Alias	BAP1; DING; BAP-1; HIP13; RING2; LUSYAM; RING1B; RING1B/RNF2
Species	Human
Gene ID (Human)	6045
Gene ID	6045
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-336 of human RING1B/RNF2 (NP_009143.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99496
Protein Weight	38kDa
Shipping	Ice bag