


Anti-Ring1B/Rnf2 Rabbit Polymab® Monoclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q99496

Item Number: CM0005892



Introduction

Rabbit PolymAb® monoclonal antibody targeting human RING1B/RNF2 (E3 ubiquitin ligase), validated for WB, IHC-P, IP, ELISA. Cross-reacts with mouse and rat. Ideal for epigenetic and chromatin biology research.

[Learn More](#)

Parameter	Value
Product Name	RING1B/RNF2 Rabbit PolymAb®
Remarks/Alias	BAP1; DING; BAP-1; HIP13; RING2; LUSYAM; RING1B
Species	Human
GeneID (Human)	6045
GeneID	6045
Immunogen	Recombinant protein of human RING1B/RNF2 (NP_009143.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99496
Protein Weight	38kDa
Shipping	Ice bag