


Anti-Riok2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bvs4

Item Number: CM0018448



Introduction

Rabbit polyclonal antibody against human RIOK2, suitable for Western blot and ELISA. Target involved in ribosome biogenesis and cell cycle regulation. Ideal for IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	RIOK2 Rabbit pAb
Remarks/Alias	RIO2; RIOK2
Species	Human
GeneID (Human)	55781
GeneID	55781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 253-552 of human RIOK2 (NP_060813.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9BVS4
Protein Weight	63kDa
Shipping	Ice bag