

Anti-Rip3/Ripk3 Polyclonal Antibody For Ihc-P And Elisa - Q9Y572

Item Number: CM0022111



Introduction

Rabbit polyclonal antibody targeting human RIP3 (RIPK3). Validated for IHC-P and ELISA. Ideal for studying necroptosis, apoptosis, and viral host interactions. Use this research-grade antibody to explore programmed cell death pathways and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	RIP3 Rabbit pAb
Remarks/Alias	RIP3; RIPK3
Species	Human
Gene ID	11035
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 410-518 of human RIP3 (NP_006862.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q9Y572
Protein Weight	57kDa
Shipping	Ice bag