


Anti-Ripk1/Rip Rabbit Monoclonal Antibody (Ko Validated) For Wb, Ip, Elisa - Q13546

Item Number: CM0004620



Introduction

KO-validated rabbit monoclonal antibody against human RIPK1/RIP. Suitable for Western blot, immunoprecipitation, and ELISA. Recombinant immunogen spans aa 465-671; detects ~76/71 kDa isoforms. Key regulator of apoptosis, necroptosis, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] RIPK1/RIP Rabbit mAb
Remarks/Alias	RIP; RIP1; AIEFL; IMD57; RIP-1
Species	Human
Gene ID (Human)	8737
Gene ID	8737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 465-671 of human RIPK1/RIP (NP_003795.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	Q13546
Protein Weight	76 kDa/71 kDa
Shipping	Ice bag