


Anti-Ripk1/Rip Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q13546

Item Number: CM0030178



Introduction

Rabbit polyclonal antibody against human RIPK1/RIP (RIP1). Suitable for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with human, mouse, rat. Research use for apoptosis, necroptosis, and inflammatory signaling.

[Learn More](#)

Parameter	Value
Product Name	RIPK1/RIP Rabbit pAb
Remarks/Alias	RIP; RIP1; AIEFL; IMD57; RIP-1; RIPK1/RIP
Species	Human
GeneID (human)	8737
GeneID	8737
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human RIPK1/RIP (NP_003795.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13546
Protein Weight	76kDa
Shipping	Ice bag