


Anti-Ripk2 Polyclonal Antibody For Wb, Elisa - 043353

Item Number: CM0019743



Introduction

High-quality Rabbit Polyclonal RIPK2 Antibody validated for WB and ELISA. Detects human, mouse, and rat RIPK2, a key innate immunity kinase involved in NOD1/2 signaling and NF- κ B activation. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	RIPK2 Rabbit pAb
Remarks/Alias	CCK; RICK; RIP2; CARD3; GIG30; CARDIAK; RIPK2
Species	Human
GeneID (Human)	8767
GeneID	8767
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-540 of human RIPK2 (NP_003812.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43353
Protein Weight	61kDa
Shipping	Ice bag