

Anti-Ripk4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P57078

Item Number: CM0028295

Introduction

CamelBio offers Anti-RIPK4 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human and mouse RIPK4. Ideal for keratinocyte differentiation, skin homeostasis, and NF-κB activation studies.

[Learn More](#)

Parameter	Value
Product Name	RIPK4 Rabbit pAb
Remarks/Alias	DIK, PKK, PPS2, RIP4, ANKK2, NKRD3, ANKRD3, CHANDS, RIPK4
Species	Human
GeneID (Human)	54101
GeneID	54101
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 240-520 of human RIPK4 (NP_065690.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P57078
Protein Weight	92kDa
Shipping	Ice bag