


Anti-Rmi2 Rabbit Polyclonal Antibody For Wb, Elisa - Q96E14

Item Number: CM0028730



Introduction

Rabbit polyclonal antibody targeting human RMI2 (BLAP18), validated for WB and ELISA. Essential for genomic stability and homologous recombination; suitable for DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	RMI2 Rabbit pAb
Remarks/Alias	BLAP18; C16orf75; RMI2
Species	Human
Gene ID	116028 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-147 of human RMI2 (NP_689521.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96E14
Protein Weight	16kDa
Shipping	Ice bag