

Anti-Rmnd1 Polyclonal Antibody For Wb, Elisa - Q9Nws8

Item Number: CM0023377

Introduction

Rabbit polyclonal antibody against human RMND1, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for mitochondrial translation and ribosome assembly research. Ideal for diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	RMND1 Rabbit pAb
Remarks/Alias	RMD1; C6orf96; COXPD11; bA351K16; bA351K16.3; RMND1
Species	Human
GeneID (Human)	55005
GeneID	55005
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 13-200 of human RMND1 (NP_060379.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NWS8
Protein Weight	52kDa
Shipping	Ice bag