

Anti-Rmnd5A Rabbit Polyclonal Antibody For Wb, Elisa - Q9H871

Item Number: CM0013733



Introduction

Rabbit polyclonal antibody against human RMND5A, a core component of the CTLH E3 ubiquitin ligase complex. Validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Supports ubiquitination pathway research.

[Learn More](#)

Parameter	Value
Product Name	RMND5A Rabbit pAb
Remarks/Alias	CTLH; GID2; RMD5; GID2A; p44CTLH; RMND5A
Species	Human
Gene ID (Human)	64795
Gene ID	64795
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 200-270 of human RMND5A (NP_073617.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H871
Protein Weight	44kDa
Shipping	Ice bag