


Anti-Rnase1 Polyclonal Antibody For Wb, Elisa - P07998

Item Number: CM0026979



Introduction

High-quality rabbit polyclonal antibody against human RNASE1 (ribonuclease A), validated for WB and ELISA, cross-reacts with mouse. Ideal for studying RNA metabolism and pancreatic function research.

[Learn More](#)

Parameter	Value
Product Name	RNASE1 Rabbit pAb
Remarks/Alias	RAC1, RIB1, RNS1, RNASE1
Species	Human
GeneID (Human)	6035
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-156 of human RNASE1 (NP_937875.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P07998
Protein Weight	18kDa
Shipping	Ice bag