

Anti-Rnase2 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P10153

Item Number: CM0030062



Introduction

Rabbit polyclonal antibody against human RNASE2 (EDN/RNS2, SWISS P10153). Validated for Western blot, immunofluorescence, and ELISA with cross-reactivity to human and mouse. Ideal for studying eosinophil function, lysosomal activity, and ribonuclease biology.

[Learn More](#)

Parameter	Specification
Product Name	RNASE2 Rabbit pAb
Remarks/Alias	EDN; RAF3; RNS2; RNASE2
Species	Human
GeneID (Human)	6036
GeneID	6036
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-161 of human RNASE2 (NP_002925.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P10153
Protein Weight	18kDa
Shipping	Ice bag