


Anti-Rnase3 Rabbit Polyclonal Antibody For Ihc-P And Elisa - P12724

Item Number: CM0023953



Introduction

RNASE3 Rabbit pAb is a polyclonal antibody targeting human eosinophil cationic protein (ECP). Validated for IHC-P and ELISA, it is suitable for innate immunity and inflammation research. - P12724

[Learn More](#)

Parameter	Value
Product Name	RNASE3 Rabbit pAb
Remarks/Alias	ECP; RAF1; RNS3; RNASE3
Species	Human
GeneID (human)	6037
GeneID	6037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-160 of human RNASE3 (NP_002926.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P12724
Protein Weight	18kDa
Shipping	Ice bag