


Anti-Rnase3 Rabbit Polyclonal Antibody For Wb, Elisa - P12724

Item Number: CM0017823



Introduction

CamelBio supplies anti-RNASE3 rabbit polyclonal antibody for WB and ELISA, targeting human eosinophil cationic protein (ECP) involved in host defense and inflammation. Cross-reacts with rat.

[Learn More](#)

Parameter	Value
Product Name	RNASE3 Rabbit pAb
Remarks/Alias	ECP; RAF1; RNS3; RNASE3
Species	Human
GeneID (Human)	6037
GeneID	6037
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 100-160 of human RNASE3 (NP_002926.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P12724
Protein Weight	18kDa
Shipping	Ice bag