


Anti-Rnase6 Polyclonal Antibody For Wb, If, Ihc, Elisa - Q93091

Item Number: CM0013931



Introduction

RNASE6 Rabbit pAb is a polyclonal antibody validated for WB, IF-P, IHC-P, and ELISA with cross-reactivity to human, mouse, and rat. Ideal for studying RNase K6 antimicrobial activity and urinary tract immunity.

[Learn More](#)

Parameter	Value
Product Name	RNASE6 Rabbit pAb
Remarks/Alias	RAD1; RNS6; RNasek6; RNASE6
Species	Human
Gene ID (Human)	6039
Gene ID	6039
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human RNASE6 (NP_005606.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q93091
Protein Weight	17kDa
Shipping	Ice bag