


Anti-Rnaseh1 Monoclonal Antibody For Wb, Elisa - O60930

Item Number: CM0006303



Introduction

Mouse monoclonal antibody targeting human RNASEH1, immunogen aa 140-286, suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for RNA-DNA hybrid and transcription studies.

[Learn More](#)

Parameter	Value
Product Name	RNASEH1 Mouse mAb
Remarks/Alias	RNH1; H1RNA; PEOB2; RNASEH1
Species	Human
GeneID	246243
Immunogen	Recombinant fusion protein containing a.a. 140-286 of human RNASEH1 (O60930)
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60930
Protein Weight	32kDa
Shipping	Ice bag