


Anti-Rnaseh1 Rabbit Polyclonal Antibody For Wb And Elisa - O60930

Item Number: CM0029118



Introduction

High-quality anti-RNASEH1 rabbit polyclonal antibody for Western blot and ELISA. Detects human, mouse, and rat RNASEH1 (32kDa). Ideal for studying RNA-DNA hybrid dynamics, transcription termination, and R-loop biology.

[Learn More](#)

Parameter	Value
Product Name	RNASEH1 Rabbit pAb
Remarks/Alias	RNH1; H1RNA; PEOB2; RNASEH1
Species	Human
Gene ID (Human)	246243
Gene ID	246243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-286 of human RNASEH1 (NP_002927.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60930
Protein Weight	32kDa
Shipping	Ice bag