


Anti-Rnaseh2A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075792

Item Number: CM0013214



Introduction

High-quality rabbit polyclonal antibody against RNASEH2A, validated for WB, IF/ICC, ELISA. Reacts with human and rat. Ideal for DNA replication and Aicardi-Goutieres syndrome research.

[Learn More](#)

Parameter	Value
Product Name	RNASEH2A Rabbit pAb
Remarks/Alias	AGS4; JUNB; RNHL; RNHIA; THSD8; RNASEHI; RNASEH2A
Species	Human
GeneID	10535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-299 of human RNASEH2A (NP_006388.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	075792
Protein Weight	33kDa
Shipping	Ice bag