


Anti-Rnaseh2C Rabbit Polyclonal Antibody For Wb And Elisa - Q8Tdp1

Item Number: CM0017620



Introduction

Rabbit polyclonal antibody targeting human RNASEH2C (AGS3/AYP1), validated for Western blot and ELISA. Ideal for studying RNA:DNA hybrid processing and DNA replication mechanisms.

[Learn More](#)

Parameter	Value
Product Name	RNASEH2C Rabbit pAb
Remarks/Alias	AGS3; AYP1; RNASEH2C
Species	Human
Gene ID (Human)	84153
Gene ID	84153
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-164 of human RNASEH2C (NP_115569.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TDP1
Protein Weight	18kDa
Shipping	Ice bag