


Anti-Rnd1 Rabbit Polyclonal Antibody For Wb - Q92730

Item Number: CM0013261



Introduction

RND1 Rabbit pAb is a polyclonal antibody targeting human RND1, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying actin cytoskeleton remodeling and neuronal process formation.

[Learn More](#)

Parameter	Description
Product Name	RND1 Rabbit pAb
Remarks/Alias	ARHS; RHO6; RHOS; RND1
Species	Human
GeneID (Human)	27289
GeneID	27289
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 133-232 of human RND1 (NP_055285.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q92730
Protein Weight	26kDa
Shipping	Ice bag