

Anti-Rnf10 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8N5U6

Item Number: CM0024157

Introduction

Rabbit polyclonal antibody targeting human RNF10, an E3 ubiquitin ligase key to ribosome quality control. Validated in WB, IF/ICC, ELISA; reacts with human, mouse, rat. Detect ~90 kDa protein, ideal for studying stress responses and proteostasis.

[Learn More](#)

Parameter	Specification
Product Name	RNF10 Rabbit pAb
Remarks/Alias	RNF10; RIE2; RING finger protein 10
Species	Human, Mouse, Rat
GeneID (human)	9921
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 660-760 of human RNF10 (NP_055683.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N5U6
Protein Weight	89kDa/90kDa
Shipping	Ice bag