

Anti-Rnf10 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8N5U6

Item Number: CM0005284

Introduction

High-quality rabbit monoclonal antibody against human RNF10, validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for studying ribosome quality control and ubiquitination in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	RNF10 Rabbit mAb
Remarks/Alias	RIE2
Species	Human
GeneID (human)	9921
GeneID	9921
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 660-760 of human RNF10 (NP_055683.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N5U6
Protein Weight	90kDa
Shipping	Ice bag