

Anti-Rnf114 Polyclonal Antibody For Wb, Elisa - Q9Y508

Item Number: CM0016851

Introduction

Rabbit polyclonal antibody against human RNF114, validated for WB and ELISA. Detects human and mouse RNF114, an E3 ubiquitin ligase involved in NF-κB regulation, cell cycle control, and innate immunity. Ideal for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	RNF114 Rabbit pAb
Remarks/Alias	ZNF313; PSORS12; RNF114
Species	Human
GeneID (Human)	55905
GeneID	55905
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human RNF114 (NP_061153.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y508
Protein Weight	26kDa
Shipping	Ice bag