

Anti-Rnf115 Polyclonal Antibody For Wb, Elisa - Q9Y4L5

Item Number: CM0020433

Introduction

Rabbit polyclonal antibody against RNF115 (BCA2, ZNF364) validated for WB and ELISA; cross-reacts with rat. RNF115 is an E3 ubiquitin ligase involved in autophagy, antiviral immunity, and cell signaling.

[Learn More](#)

Parameter	Value
Product Name	RNF115 Rabbit pAb
Remarks/Alias	BCA2; ZFP364; ZNF364; RABRING7; RNF115
Species	Human
GeneID (Human)	27246
GeneID	27246
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 30-140 of human RNF115 (NP_055270.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9Y4L5
Protein Weight	34kDa
Shipping	Ice bag