


Anti-Rnf125 Polyclonal Antibody For Wb, Elisa - Q96Eq8

Item Number: CM0015317



Introduction

RNF125 Rabbit pAb is a rabbit polyclonal antibody validated for Western blot and ELISA, targeting human RNF125, an E3 ubiquitin ligase regulating antiviral immunity and T-cell activation. Suitable for research applications.

[Learn More](#)

Parameter	Value
Product Name	RNF125 Rabbit pAb
Remarks/Alias	TNORS; TRAC1; TRAC-1; RNF125
Species	Human
GeneID	54941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 41-200 of human RNF125 (NP_060301.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96EQ8
Protein Weight	26kDa
Shipping	Ice bag