


Anti-Rnf126 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bv68

Item Number: CM0021451



Introduction

Rabbit polyclonal antibody against RNF126, validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for studying E3 ubiquitin ligase functions in protein quality control and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	RNF126 Rabbit pAb
Remarks/Alias	RNF126; ring finger protein 126
Species	Human
Gene ID (Human)	55658
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-200 of human RNF126 (NP_919442.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV68
Protein Weight	33kDa/35kDa
Shipping	Ice bag