


Anti-Rnf126 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Bv68

Item Number: CM0019501



Introduction

Rabbit polyclonal antibody against human E3 ubiquitin-protein ligase RNF126, suitable for Western blot, IHC-P, IF/ICC, and ELISA. Detects 34 kDa RNF126 in human, mouse, and rat for ubiquitination and DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	RNF126 Rabbit pAb
Remarks/Alias	RNF126
Species	Human
Gene ID (Human)	55658
Gene ID	55658
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV68
Protein Weight	34kDa
Shipping	Ice bag