

Anti-Rnf14 Monoclonal Antibody For Wb, Elisa - Q9Ubs8

Item Number: CM0003601

Introduction

Anti-RNF14 rabbit monoclonal antibody validated for WB and ELISA. Detects human, mouse, and rat RNF14, a key E3 ligase in translation quality control and Wnt signaling. Ideal for mechanistic and diagnostic research.

[Learn More](#)

Specification	Details
Product Name	RNF14 Rabbit mAb
Remarks/Alias	ARA54; HFB30; TRIAD2; HRIHFB2038; RNF14
Species	Human
Gene ID (Human)	9604
Gene ID	9604
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBS8
Protein Weight	54kDa
Shipping	Ice bag