


Anti-Rnf144A Polyclonal Antibody For Wb, Elisa - P50876

Item Number: CM0020243



Introduction

Rabbit polyclonal antibody targeting human RNF144A (UIP4), validated for Western blot and ELISA. Cross-reacts with rat. Suitable for ubiquitination, DNA damage, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	RNF144A Rabbit pAb
Remarks/Alias	UIP4; hUIP4; RNF144; UBCE7IP4
Species	Human
GeneID (Human)	9781
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 52-179 of human RNF144A (NP_055561.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P50876
Protein Weight	33kDa
Shipping	Ice bag