

Anti-Rnf149 Polyclonal Antibody For If/Icc, Elisa - Q8Nc42

Item Number: CM0024536

Introduction

Rabbit polyclonal antibody targeting human RNF149 (43kDa) for IF/ICC and ELISA. Detects endogenous RNF149; cross-reactivity with rat. Ideal for E3 ligase and BRAF degradation research.

[Learn More](#)

Parameter	Specification
Product Name	RNF149 Rabbit pAb
Remarks/Alias	RNF149; DNAPTP2; ring finger protein 149
Species	Human
GeneID (Human)	284996
GeneID	284996
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human RNF149(NP_775918.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8NC42
Protein Weight	43kDa
Shipping	Ice bag