

Anti-Rnf149 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q8Nc42

Item Number: CM0015386

Introduction

Rabbit polyclonal antibody against human RNF149, validated for WB, IHC-P, IF/ICC, and ELISA. Detects RNF149 in human, mouse, and rat samples. Useful for studying ubiquitination and BRAF degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	RNF149 Rabbit pAb
Remarks/Alias	DNATP2; RNF149
Species	Human
Gene ID (Human)	284996
Gene ID	284996
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 221-400 of human RNF149 (NP_775918.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NC42
Protein Weight	43kDa
Shipping	Ice bag