

Anti-Rnf152 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N8N0

Item Number: CM0019657

Introduction

Anti-RNF152 rabbit polyclonal antibody for WB and ELISA, targeting human E3 ubiquitin-protein ligase RNF152 (UniProt Q8N8N0). Validated for research on mTORC1 signaling and apoptosis pathways.

[Learn More](#)

Parameter	Value
Product Name	RNF152 Rabbit pAb
Remarks/Alias	RNF152
Species	Human
Gene ID (Human)	220441
Gene ID	220441
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human RNF152 (NP_775828.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N8N0
Protein Weight	22kDa
Shipping	Ice bag