

Anti-Rnf166 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96A37

Item Number: CM0017194

Introduction

Rabbit polyclonal antibody specific for human RNF166, validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse. Targets E3 ubiquitin ligase involved in autophagy and innate immunity. Suitable for IVD research applications.

[Learn More](#)

Parameter	Value
Product Name	RNF166 Rabbit pAb
Remarks/Alias	RNF166
Species	Human
GeneID (Human)	115992
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-237 of human RNF166 (NP_849163.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96A37
Protein Weight	26kDa
Shipping	Ice bag