


Anti-Rnf168 Polyclonal Antibody For Wb, Elisa - Q8Iyw5

Item Number: CM0030001



Introduction

Rabbit polyclonal antibody against human E3 ubiquitin-protein ligase RNF168, validated for Western blotting and ELISA. Detects human and mouse RNF168, a key regulator of DNA damage response and histone ubiquitination.

[Learn More](#)

Parameter	Value
Product Name	RNF168 Rabbit pAb
Remarks/Alias	RIDL; hRNF168; RNF168
Species	Human
GeneID (Human)	165918
GeneID	165918
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 300-571 of human RNF168 (NP_689830.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8IYW5
Protein Weight	65kDa
Shipping	Ice bag