

Anti-Rnf180 Rabbit Polyclonal Antibody For Wb And Elisa - Q86T96

Item Number: CM0020373

Introduction

Rabbit polyclonal antibody raised against recombinant human RNF180 (aa 1-416). Recognizes human, mouse, rat RNF180. Validated for western blot and ELISA. Ideal for ubiquitination and protein degradation studies.

[Learn More](#)

Parameter	Specification
Product Name	RNF180 Rabbit pAb
Remarks/Alias	RINES; RNF180
Species	Human
Gene ID (Human)	285671
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-416 of human RNF180 (NP_848627.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86T96
Protein Weight	68kDa
Shipping	Ice bag