

Anti-Rnf181 Rabbit Polyclonal Antibody For Wb And Elisa - Q9POP0

Item Number: CM0015889

Introduction

Anti-RNF181 rabbit polyclonal antibody validated for WB and ELISA. Reacts with human, mouse, rat. Target is E3 ubiquitin-protein ligase RNF181, involved in ubiquitination and proteasome regulation. Ideal for ubiquitin pathway research.

[Learn More](#)

Parameter	Value
Product Name	RNF181 Rabbit pAb
Remarks/Alias	HSPC238; RNF181
Species	Human
Gene ID	51255
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human RNF181 (NP_057578.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9POP0
Protein Weight	18kDa
Shipping	Ice bag