

Anti-Rnf185 Rabbit Polyclonal Antibody For Wb And Elisa - Q96Gf1

Item Number: CM0024401

Introduction

Rabbit polyclonal antibody against human RNF185, suitable for WB and ELISA. Cross-reacts with mouse and rat. Useful for research on mitophagy, ER stress, ERAD, and innate antiviral immunity.

[Learn More](#)

Parameter	Value
Product Name	RNF185 Rabbit pAb
Remarks/Alias	RNF185
Species	Human
GeneID (Human)	91445
GeneID	91445
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human RNF185 (NP_689480.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96GF1
Protein Weight	20kDa
Shipping	Ice bag