


Anti-Rnf20 Rabbit Monoclonal Antibody For Wb And Elisa - Q5Vtr2

Item Number: CM0005904



Introduction

High-quality anti-RNF20 rabbit monoclonal antibody for Western blot and ELISA. Reacts with human and mouse RNF20, a histone H2B ubiquitin ligase. Ideal for epigenetic and cancer research.

[Learn More](#)

Parameter	Value
Product Name	RNF20 Rabbit mAb
Remarks/Alias	BRE1; BRE1A; hBRE1; RNF20
Species	Human
Gene ID (Human)	56254
Gene ID	56254
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1-100 of human RNF20 (Q5VTR2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5VTR2
Protein Weight	114kDa
Shipping	Ice bag