

Anti-Rnf220 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Vtb9

Item Number: CM0013694

Introduction

High-quality rabbit polyclonal antibody against RNF220 (E3 ubiquitin-protein ligase). Validated for Western blotting and ELISA. Recognizes human, mouse, and rat RNF220 (63kDa). Ideal for studying ubiquitination, Wnt signaling, neurodevelopment, and interferon responses.

[Learn More](#)

Parameter	Value
Product Name	RNF220 Rabbit pAb
Remarks/Alias	HLD23; C1orf164; RNF220
Species	Human
Gene ID (Human)	55182
Gene ID	55182
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 267-566 of human RNF220 (NP_060620.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5VTB9
Protein Weight	63kDa
Shipping	Ice bag