

Anti-Rnf31 Polyclonal Antibody For Wb, If, Elisa - Q96Ep0

Item Number: CM0016423



Introduction

Rabbit polyclonal antibody targeting human RNF31 (HOIP, ZIBRA) for Western blotting, immunofluorescence, and ELISA. Validated in human and mouse samples, recognizing the ~120 kDa E3 ligase component of the LUBAC complex. Ideal for NF- κ B signaling and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	RNF31 Rabbit pAb
Remarks/Alias	HOIP; Paul; ZIBRA; RNF31
Species	Human
GeneID (Human)	55072
GeneID	55072
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human RNF31 (NP_060469.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96EP0
Protein Weight	120kDa
Shipping	Ice bag