

Anti-Rnf31 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Ep0

Item Number: CM0026720



Introduction

High-quality rabbit polyclonal antibody against RNF31 (HOIP), validated for WB and ELISA. Reacts with human and rat. Ideal for studying NF-κB signaling and linear ubiquitination.

[Learn More](#)

Parameter	Value
Product Name	RNF31 Rabbit pAb
Remarks/Alias	HOIP; Paul; ZIBRA; RNF31
Species	Human
GeneID (Human)	55072
GeneID	55072
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 500-700 of human RNF31 (NP_060469.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q96EP0
Protein Weight	120kDa
Shipping	Ice bag