

Anti-Rnf34 Rabbit Polyclonal Antibody For Wb, If/Icc, And Elisa - Q969K3

Item Number: CM0027248



Introduction

High-quality rabbit polyclonal antibody targeting human RNF34 (CARP1/hRFI), an E3 ubiquitin ligase. Validated for WB, IF/ICC, and ELISA applications. Ideal for apoptosis, protein degradation, and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	RNF34 Rabbit pAb
Remarks/Alias	RFI, RIF, RIFF, hRFI, CARP1, CARP-1, RNF34
Species	Human
GeneID	80196
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human RNF34 (NP_079402.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q969K3
Protein Weight	42kDa
Shipping	Ice bag