

Anti-Rnf4 Rabbit Polyclonal Antibody For Wb, Elisa - P78317

Item Number: CM0013519



Introduction

Rabbit polyclonal antibody against human RNF4 (SLX5/SNURF), validated for WB and ELISA, with cross-reactivity to mouse and rat. Targets the E3 ubiquitin-protein ligase involved in sumoylation-dependent degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	RNF4 Rabbit pAb
Remarks / Alias	SLX5, SNURF, RES4-26, RNF4
Species	Human
Gene ID	6047
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human RNF4 (NP_001171938.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P78317
Protein Weight	21 kDa
Shipping	Ice bag